

Anti-A-Smooth Muscle Actin (Acta2) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P62736

Item Number: CM0013516

Introduction

High-quality rabbit polyclonal antibody against human α -Smooth Muscle Actin (ACTA2, P62736). Validated for WB, IHC-P, IF-P, ELISA. Cross-reacts with mouse, rat. Ideal for cytoskeleton and smooth muscle research.

[Learn More](#)

Parameter	Value
Product Name	α -Smooth Muscle Actin (ACTA2) Rabbit pAb
Remarks/Alias	ACTSA; α -Smooth Muscle Actin (ACTA2)
Species	Human
Gene ID (Human)	59
Gene ID	59
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human α -Smooth Muscle Actin (ACTA2) (NP_001604.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62736
Protein Weight	42kDa
Shipping	Ice bag