


Anti-A-Synuclein Monoclonal Antibody For Wb - P37840

Item Number: CM0006519



Introduction

Rabbit monoclonal antibody against human α -Synuclein (PD1, NACP, PARK1/4). Validated for WB, IF/ICC, IF-P, ELISA; cross-reacts with mouse and rat. Essential for synaptic function and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	α -Synuclein Rabbit mAb
Remarks/Alias	PD1; NACP; PARK1; PARK4; α -Synuclein
Species	Human
Gene ID (Human)	6622
Gene ID	6622
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 59-140 of human α -Synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag