


Anti-A-Synuclein Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - P37840

Item Number: CM0006570



Introduction

Recombinant α -Synuclein rabbit monoclonal antibody (P37840) validated for WB, IF, IHC, and ELISA. Reacts with human, mouse, rat. 14 kDa. Suitable for synaptic function, neurodegeneration, and Parkinson's disease research.

[Learn More](#)

Parameter	Details
Product Name	α -Synuclein Rabbit mAb
Remarks/Alias	PD1, NACP, PARK1, PARK4, α -Synuclein
Species	Human
Gene ID	6622 (Human)
Immunogen	Recombinant protein corresponding to amino acids 1-140 of human α -Synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14 kDa
Shipping	Ice bag