


Anti-A-Synuclein Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Elisa - P37840

Item Number: CM0007447



Introduction

High-quality rabbit monoclonal antibody against human α -synuclein (P37840), validated for WB, IP, IF/ICC, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for neurodegeneration research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	α -Synuclein Rabbit mAb
Remarks/Alias	PD1, NACP, PARK1, PARK4, α -Synuclein
Species	Human
GeneID (Human)	6622
GeneID	6622
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human α -Synuclein (NP_000336.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37840
Protein Weight	14kDa
Shipping	Ice bag