

Anti-A-Tubulin Rabbit Monoclonal Antibody For Wb - P68366

Item Number: CM0011002

Introduction

High-quality anti- α -tubulin rabbit monoclonal antibody validates for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat α -tubulin (~50 kDa). Ideal for microtubule and cytoskeleton research studies.

[Learn More](#)

Parameter	Value
Product Name	α -Tubulin Rabbit mAb
Remarks/Alias	ALS22; TUBA1; H2-ALPHA; α -Tubulin
Species	Human
Gene ID (Human)	7277
Gene ID	7277
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 349-448 of human α -Tubulin (P68366).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68366
Protein Weight	50 kDa
Shipping	Ice bag