

Anti-A1Bg Polyclonal Antibody For Wb, Elisa - P04217

Item Number: CM0014304

Introduction

A1BG Rabbit pAb (Polyclonal Antibody) validated for WB, ELISA with Human, Mouse, Rat cross-reactivity. Recombinant immunogen, 54kDa. Ideal for plasma protein research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	A1BG Rabbit pAb
Remarks/Alias	A1B; ABG; GAB; HYST2477; A1BG
Species	Human
Gene ID (Human)	1
Gene ID	1
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 236-495 of human A1BG (NP_570602.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04217
Protein Weight	54kDa
Shipping	Ice bag