

Anti-A1Cf (Apobec1 Complementation Factor) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nq94

Item Number: CM0014005



Introduction

A1CF rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat APOBEC1 complementation factor, an essential component of the APOB mRNA editing complex.

[Learn More](#)

Parameter	Value
Product Name	A1CF Rabbit pAb
Remarks/Alias	ACF; ASP; ACF64; ACF65; APOBEC1CF; A1CF
Species	Human
Gene ID (Human)	29974
Gene ID	29974
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-120 of human A1CF (NP_055391.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ94
Protein Weight	65kDa
Shipping	Ice bag