


Anti-A2ML1 Rabbit Polyclonal Antibody For Wb, Elisa - A8K2U0

Item Number: CM0024630



Introduction

High-quality rabbit polyclonal antibody targeting human A2ML1, validated for Western blot and ELISA. Shows cross-reactivity with rat. Ideal for investigating protease inhibition mechanisms and epidermal biology, including desquamation processes.

[Learn More](#)

Parameter	Value
Product Name	A2ML1 Rabbit pAb
Remarks/Alias	OMS; p170; CPAMD9; A2ML1
Species	Human
Gene ID (Human)	144568
Gene ID	144568
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human A2ML1 (NP_653271.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	A8K2U0
Protein Weight	161kDa
Shipping	Ice bag