


Anti-Aaas Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nrg9

Item Number: CM0029194



Introduction

AAAS rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Reacts with human Aladin (AAAS) protein (Gene ID: 8086, ~60 kDa). Ideal for research on nuclear pore complex function, nervous system development, and mitosis.

[Learn More](#)

Parameter	Value
Product Name	AAAS Rabbit pAb
Remarks/Alias	AAA, AAASb, GL003, ALADIN, ADRACALA, ADRACALIN, AAAS
Species	Human
Gene ID	8086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 287-546 of human AAAS (NP_056480.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NRG9
Protein Weight	60kDa
Shipping	Ice bag