

Anti-Aadac Polyclonal Antibody For Wb, Elisa - P22760

Item Number: CM0015677

Introduction

Rabbit polyclonal antibody targeting human AADAC (arylacetamide deacetylase), validated for Western blot and ELISA. Cross-reacts with mouse. Recombinant immunogen (aa 170-399). Useful for studying triglyceride lipase activity, xenobiotic metabolism, and liver-related pathways.

[Learn More](#)

Parameter	Value
Product Name	AADAC Rabbit pAb
Remarks/Alias	DAC; CES5A1; AADAC
Species	Human
GeneID (Human)	13
GeneID	13
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-399 of human AADAC (NP_001077.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22760
Protein Weight	46kDa
Shipping	Ice bag