

Anti-Aadat Rabbit Polyclonal Antibody For Wb, Elisa - Q8N5Z0

Item Number: CM0018322



Introduction

AADAT Rabbit pAb is a polyclonal antibody raised in rabbit against human AADAT. Suitable for WB and ELISA, it cross-reacts with rat and aids research on aminotransferases and kynurenine metabolism.

[Learn More](#)

Parameter	Value
Product Name	AADAT Rabbit pAb
Remarks/Alias	KAT2; KATII; KYAT2; AADAT
Species	Human
GeneID (Human)	51166
GeneID	51166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-429 of human AADAT (NP_001273611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q8N5Z0
Protein Weight	47kDa
Shipping	Ice bag