

Anti-Aadat Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N5Z0

Item Number: CM0013379



Introduction

AADAT Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA applications. It detects AADAT (KAT2/KYAT2) in human, mouse, and rat samples, ideal for studying mitochondrial aminotransferase functions.

[Learn More](#)

Parameter	Value
Product Name	AADAT Rabbit pAb
Remarks/Alias	KAT2; KATII; KYAT2; AADAT
Species	Human
Gene ID (Human)	51166
Gene ID	51166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-429 of human AADAT (NP_001273611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5Z0
Protein Weight	47kDa
Shipping	Ice bag