


Anti-Aak1 Rabbit Polyclonal Antibody For Wb, Elisa - Q2M2I8

Item Number: CM0024199



Introduction

Rabbit polyclonal antibody against human AAK1, validated for Western blotting and ELISA. Targets AP2-associated protein kinase 1, a regulator of clathrin-mediated endocytosis. Ideal for research on endocytosis and viral entry.

[Learn More](#)

Parameter	Value
Product Name	AAK1 Rabbit pAb
Species	Human
Gene ID	22848
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 450-550 of human AAK1 (NP_055726.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q2M2I8
Protein Weight	103 kDa
Shipping	Ice bag