

## Anti-Aamp Monoclonal Antibody For Wb, If/Icc, Elisa - Q13685

Item Number: CM0002613

### Introduction

Anti-AAMP rabbit monoclonal antibody validated for Western blot (WB), immunofluorescence/immunocytochemistry (IF/ICC), and ELISA; cross-reacts with human, mouse, and rat AAMP. Ideal for angiogenesis and cell migration studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | AAMP Rabbit mAb                                                                                                     |
| Remarks/Alias    | AAMP                                                                                                                |
| Species          | Human                                                                                                               |
| GeneID (Human)   | 14                                                                                                                  |
| GeneID           | 14                                                                                                                  |
| Immunogen        | Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human AAMP (Q13685). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q13685                                                                                                              |
| Protein Weight   | 47kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |