


Anti-Aanat Rabbit Polyclonal Antibody For Wb And Elisa - Q16613

Item Number: CM0014228



Introduction

Rabbit polyclonal antibody targeting human serotonin N-acetyltransferase (AANAT) for Western blotting and ELISA. Recombinant protein immunogen (aa 1-207). Useful for melatonin biosynthesis and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	AANAT Rabbit pAb
Remarks/Alias	DSPS; SNAT; AANAT
Species	Human
Gene ID (Human)	15
Gene ID	15
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human AANAT (NP_001079.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16613
Protein Weight	23kDa
Shipping	Ice bag