

Anti-Aar2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y312

Item Number: CM0019338



Introduction

High-quality anti-AAR2 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human and mouse AAR2, a U5 snRNP component essential for spliceosome assembly and pre-mRNA splicing.

[Learn More](#)

Parameter	Value
Product Name	AAR2 Rabbit pAb
Remarks/Alias	CGI-23; C20orf4; AAR2
Species	Human
GeneID (Human)	25980
GeneID	25980
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human AAR2 (NP_056326.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y312
Protein Weight	43kDa
Shipping	Ice bag