

Anti-Aars2 Polyclonal Antibody For Wb, Elisa - Q5Jtz9

Item Number: CM0027023



Introduction

Rabbit polyclonal antibody targeting human AARS2 (Alanyl-tRNA synthetase, mitochondrial), suitable for WB and ELISA. Ideal for studying mitochondrial protein biosynthesis, editing, and protein lactylation. - Q5JTZ9

[Learn More](#)

Parameter	Value
Product Name	AARS2 Rabbit pAb
Remarks/Alias	AARSL; LKENP; COXPD8; MTALARS; MT-ALARS; AARS2
Species	Human
Gene ID (Human)	57505
Gene ID	57505
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 841-985 of human AARS2 (NP_065796.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5JTZ9
Protein Weight	107kDa
Shipping	Ice bag