


Anti-Aars2 Rabbit Polyclonal Antibody For Wb And Elisa - Q5Jtz9

Item Number: CM0026210



Introduction

Validated AARS2 rabbit polyclonal antibody for WB and ELISA applications, detecting human and mouse mitochondrial alanine-tRNA ligase. Involved in protein translation and cGAS/STING signaling inhibition.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] AARS2 Rabbit pAb
Remarks/Alias	AARSL; LKENP; COXPD8; MTALARS; MT-ALARS
Species	Human
Gene ID	57505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 491-620 of human AARS2 (NP_065796.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5J TZ9
Protein Weight	107kDa
Shipping	Ice bag