

Anti-Aars/Aars1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P49588

Item Number: CM0015108



Introduction

Rabbit polyclonal antibody against human alanyl-tRNA synthetase (AARS/AARS1). Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Targets nuclear and cytoplasmic forms; useful for protein lactylation and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	AARS Rabbit pAb
Remarks/Alias	AARS; TTD8; CMT2N; DEE29; HDLS2; EIEE29
Species	Human
Gene ID	16
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49588
Protein Weight	107 kDa
Shipping	Ice bag