


Anti-Aatf Polyclonal Antibody For Wb, Elisa - Q9Ny61

Item Number: CM0025540



Introduction

Rabbit polyclonal antibody targeting AATF (Che-1/DED), validated for Western blotting and ELISA. Reacts with human and mouse. AATF is a nucleolar protein involved in ribosome biogenesis and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	AATF Rabbit pAb
Remarks/Alias	DED; BFR2; CHE1; CHE-1; AATF
Species	Human
Gene ID (Human)	26574
Gene ID	26574
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 7-100 of human AATF (NP_036270.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NY61
Protein Weight	63kDa
Shipping	Ice bag