

Anti-Abca1 Polyclonal Antibody For Wb, If/Icc, Elisa - 095477

Item Number: CM0022051



Introduction

ABCA1 Rabbit pAb detects human ABCA1 (TGD, ABC1, CERP) with cross-reactivity to mouse. Suitable for WB, IF/ICC, ELISA. Ideal for cholesterol efflux and HDL formation studies.

[Learn More](#)

Parameter	Value
Product Name	ABCA1 Rabbit pAb
Remarks/Alias	TGD; ABC1; CERP; ABC-1; HDLDT1; HPALP1; HDLCQTL13; ABCA1
Species	Human
Gene ID (Human)	19
Gene ID	19
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1170-1350 of human ABCA1 (NP_005493.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	095477
Protein Weight	254 kDa
Shipping	Ice bag