

Anti-Abca1 Rabbit Monoclonal Antibody For Wb, Elisa - 095477

Item Number: CM0012628



Introduction

ABCA1 Rabbit Monoclonal Antibody validated for Western blotting (WB) and ELISA, with cross-reactivity to human, mouse, and rat. Targets phospholipid-transporting ATPase ABCA1, a key regulator of cholesterol efflux and HDL formation, making it suitable for lipid metabolism and atherosclerosis research.

[Learn More](#)

Parameter	Value
Product Name	ABCA1 Rabbit mAb
Remarks/Alias	TGD; ABC1; CERP; ABC-1; HDLDT1; HPALP1; HDLCQTL13
Species	Human
GeneID (Human)	19
GeneID	19
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1170-1350 of human ABCA1 (NP_005493.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095477
Protein Weight	254kDa
Shipping	Ice bag