


Anti-Abca1 Rabbit Monoclonal Antibody For Wb, Elisa - O95477

Item Number: CM0007168



Introduction

High-quality rabbit monoclonal antibody against human ABCA1, validated for WB and ELISA. Detects endogenous ABCA1 (254 kDa) in human, mouse, and rat samples. Ideal for lipid metabolism and HDL research.

[Learn More](#)

Parameter	Value
Product Name	ABCA1 Rabbit mAb
Remarks/Alias	TGD; ABC1; CERP; ABC-1; HDLDT1; HPALP1; HDLCQTL13; ABCA1
Species	Human
Gene ID (Human)	19
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1170-1350 of human ABCA1 (NP_005493.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95477
Protein Weight	254kDa
Shipping	Ice bag