

Anti-Abca1 Rabbit Polyclonal Antibody - 095477

Item Number: CM0025277



Introduction

High-quality rabbit polyclonal antibody against human ABCA1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for cholesterol efflux and HDL research.

[Learn More](#)

Parameter	Value
Product Name	ABCA1 Rabbit pAb
Remarks/Alias	TGD, ABC1, CERP, ABC-1, HDLDT1, HPALP1, HDLCQTL13, ABCA1
Species	Human
GeneID (Human)	19
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1870-2120 of human ABCA1 (NP_005493.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95477
Protein Weight	254 kDa
Shipping	Ice bag