

Anti-Abca12 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q86Uk0

Item Number: CM0021634

Introduction

High-quality ABCA12 Rabbit Polyclonal Antibody validated for immunofluorescence (IF/ICC) and ELISA, with mouse cross-reactivity. Targets the glucosylceramide transporter ABCA12, essential for epidermal lipid lamellae formation and skin barrier function. Ideal for research in dermatology, lipid transport, and ichthyosis.

[Learn More](#)

Parameter	Value
Product Name	ABCA12 Rabbit pAb
Remarks/Alias	LI2; ICR2B; ARCI4A; ARCI4B; ABCA12
Species	Human
Gene ID (Human)	26154
Gene ID	26154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2456-2595 of human ABCA12 (NP_775099.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q86UK0
Protein Weight	293kDa
Shipping	Ice bag