

Anti-Abca2 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Bzc7

Item Number: CM0024636

Introduction

Rabbit polyclonal antibody against human ABCA2, validated for WB, IF/ICC, IF-P, ELISA. Reacts with human, mouse, rat. ABCA2 is a lipid transporter involved in cholesterol homeostasis and myelin formation.

[Learn More](#)

Parameter	Value
Product Name	ABCA2 Rabbit pAb
Remarks/Alias	ABC2; IDPOGSA; ABCA2
Species	Human
Gene ID (Human)	20
Gene ID	20
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-350 of human ABCA2 (NP_001597.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZC7
Protein Weight	270kDa
Shipping	Ice bag