

Anti-Abca4 Rabbit Polyclonal Antibody For Wb, Elisa - P78363

Item Number: CM0012846



Introduction

ABCA4 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. It targets the retinal-specific ABC transporter ABCA4, essential for retinoid clearance in photoreceptors.

[Learn More](#)

Parameter	Value
Product Name	ABCA4 Rabbit pAb
Remarks/Alias	FFM; RMP; ABCR; RP19; STGD; ABC10; ARMD2; CORD3; STGD1; ABCA4
Species	Human
GeneID (Human)	24
GeneID	24
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 110-300 of human ABCA4 (NP_000341.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78363
Protein Weight	256kDa
Shipping	Ice bag