


Anti-Abca7 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Izy2

Item Number: CM0023778



Introduction

Rabbit polyclonal antibody targeting human ABCA7, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Phospholipid transporter implicated in Alzheimer's disease and lipid homeostasis.

[Learn More](#)

Parameter	Value
Product Name	ABCA7 Rabbit pAb
Remarks/Alias	AD9; ABCX; ABCA-SSN; ABCA7
Species	Human
Gene ID (Human)	10347
Gene ID	10347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2010-2146 of human ABCA7 (NP_061985.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZY2
Protein Weight	234kDa
Shipping	Ice bag