

Anti-Abca8A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8K442

Item Number: CM0021678

Introduction

Rabbit polyclonal antibody targeting mouse ABCA8A, suitable for Western blot, immunohistochemistry, and ELISA. Cross-reacts with human and rat. Useful for cholesterol efflux and HDL regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ABCA8A Rabbit pAb
Species	Mouse
Gene ID (Mouse)	217258
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 690-828 of mouse ABCA8A (NP_694785.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8K442
Protein Weight	184kDa
Shipping	Ice bag