


Anti-Abcb11 Rabbit Polyclonal Antibody For If/Icc And Elisa - O95342

Item Number: CM0018732



Introduction

High-quality rabbit polyclonal antibody against human ABCB11 (BSEP), validated for IF/ICC, ELISA. Detects mouse and rat. Ideal for bile salt export pump research, hepatic transport studies, and cholestasis investigation.

[Learn More](#)

Parameter	Value
Product Name	ABCB11 Rabbit pAb
Remarks/Alias	BSEP; PGY4; SPGP; ABC16; BRIC2; PFIC2; PFIC-2; ABCB11
Species	Human
GeneID (Human)	8647
GeneID	8647
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1050-1150 of human ABCB11 (NP_003733.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95342
Protein Weight	146kDa
Shipping	Ice bag