


Anti-Abcb11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O95342

Item Number: CM0027707



Introduction

High-quality rabbit polyclonal antibody against ABCB11 (BSEP), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for studies of bile salt transport and intrahepatic cholestasis.

[Learn More](#)

Parameter	Value
Product Name	ABCB11 Rabbit pAb
Remarks/Alias	BSEP; PGY4; SPGP; ABC16; BRIC2; PFIC2; PFIC-2; ABCB11
Species	Human
Gene ID (Human)	8647
Gene ID	8647
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1050-1150 of human ABCB11 (NP_003733.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95342
Protein Weight	146kDa
Shipping	Ice bag