

Anti-Abcb4 Polyclonal Antibody For Wb, Elisa - P21439

Item Number: CM0030829



Introduction

High-quality anti-ABCB4 polyclonal antibody validated for WB and ELISA applications. Reacts with human and mouse ABCB4, a phospholipid floppase crucial for biliary lipid secretion. Ideal for hepatobiliary disease research.

[Learn More](#)

Parameter	Value
Product Name	ABCB4 Rabbit pAb
Remarks/Alias	GBD1; ICP3; MDR2; MDR3; PGY3; ABC21; MDR2/3; PFIC-3; ABCB4
Species	Human
GeneID (Human)	5244
GeneID	5244
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human ABCB4 (NP_061337.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21439
Protein Weight	142kDa
Shipping	Ice bag