


Anti-Abcb7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75027

Item Number: CM0015236



Introduction

High-quality rabbit polyclonal antibody against human ABCB7, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial iron-sulfur cluster transport and heme biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	ABCB7 Rabbit pAb
Remarks/Alias	ABCB7; ASAT; Atm1p; EST140535; ABCB7
Species	Human
Gene ID	22
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 503-753 of human ABCB7 (NP_004290.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75027
Protein Weight	83kDa
Shipping	Ice bag