


Anti-Abcb8 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nut2

Item Number: CM0030529



Introduction

Rabbit polyclonal antibody against human ABCB8, validated for WB, IHC-P, IF/ICC, ELISA. Target the mitochondrial potassium channel ATP-binding subunit involved in iron transport and cardiac function. Cross-reacts with mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	ABCB8 Rabbit pAb
Remarks/Alias	MABC1; M-ABC1; MITOSUR; EST328128; ABCB8
Species	Human
Gene ID (Human)	11194
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 469-718 of human ABCB8 (NP_009119.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NUT2
Protein Weight	80kDa
Shipping	Ice bag