


Anti-Abcc5 Rabbit Polyclonal Antibody For Wb, Elisa - O15440

Item Number: CM0025210



Introduction

ABCC5 Rabbit pAb validated for WB, ELISA. Detects human, mouse, rat ABCC5 (MRP5). Recombinant immunogen (aa 1-179). For transport and drug resistance research.

[Learn More](#)

Parameter	Value
Product Name	ABCC5 Rabbit pAb
Remarks/Alias	MRP5; SMRP; ABC33; MOATC; MOAT-C; pABC11; EST277145; ABCC5
Species	Human
GeneID (Human)	10057
GeneID	10057
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human ABCC5 (NP_005679.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15440
Protein Weight	161kDa
Shipping	Ice bag