

Anti-Abcc5/Mrp5 Rabbit Polyclonal Antibody For If-P And Elisa - O15440

Item Number: CM0025839



Introduction

Rabbit polyclonal antibody targeting human ABCC5/MRP5, validated for IF-P and ELISA applications. MRP5 is an ATP-dependent transporter implicated in drug resistance and metabolite efflux. Suitable for research in oncology and pharmacology.

[Learn More](#)

Parameter	Specification
Product Name	ABCC5/MRP5 Rabbit pAb
Remarks/Alias	MRP5; SMRP; ABC33; MOATC; MOAT-C; pABC11; EST277145
Species	Human
GeneID	10057
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1101-1200 of human ABCC5/MRP5 (NP_005679.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Human
SWISS	O15440
Protein Weight	161kDa
Shipping	Ice bag