

Anti-Abcc6 Rabbit Polyclonal Antibody For Wb, Elisa - O95255

Item Number: CM0027206



Introduction

Rabbit polyclonal antibody targeting human ABCC6, suitable for WB and ELISA. Shows cross-reactivity with mouse and rat. Ideal for studying ABC transporters, drug resistance, and pseudoxanthoma elasticum.

[Learn More](#)

Parameter	Value
Product Name	ABCC6 Rabbit pAb
Remarks/Alias	ARA; PXE; MLP1; MRP6; PXE1; URG7; ABC34; GAC12; MOATE; MOAT-E; EST349056; ABCC6
Species	Human
GeneID (Human)	368
GeneID	368
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 830-940 of human ABCC6 (NP_001162.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95255
Protein Weight	165kDa
Shipping	Ice bag