


Anti-Abcd1 / Ald Rabbit Monoclonal Antibody For Wb, Elisa - P33897

Item Number: CM0007568



Introduction

Recombinant rabbit monoclonal antibody against ABCD1/ALDP (amino acids 650-745). Suitable for Western blot and ELISA; cross-reacts with human and mouse. ABCD1 is an ATP-dependent VLCFA transporter critical for peroxisomal fatty acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ABCD1 / ALD Rabbit mAb
Remarks/Alias	ALD; AMN; ALDP; ABC42; ABCD1 / ALD
Species	Human
Gene ID (Human)	215
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 650-745 of human ABCD1 / ALD (NP_000024.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P33897
Protein Weight	83kDa
Shipping	Ice bag