


Anti-Abcd1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P33897

Item Number: CM0013027



Introduction

Rabbit polyclonal antibody targeting human ABCD1 (ALDP), validated for WB, IF/ICC, and ELISA. Detects human and mouse ABCD1; suitable for peroxisomal transport and adrenoleukodystrophy research.

[Learn More](#)

Parameter	Value
Product Name	ABCD1 Rabbit pAb
Remarks/Alias	ALD; AMN; ALDP; ABC42; ABCD1
Species	Human
Gene ID (Human)	215
Gene ID	215
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 340-500 of human ABCD1 (NP_000024.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P33897
Protein Weight	83kDa
Shipping	Ice bag