


Anti-Abce1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P61221

Item Number: CM0006184



Introduction

Anti-ABCE1 rabbit monoclonal antibody for WB, IF/ICC, ELISA. Cross-reacts with human, mouse. Targets ribosome recycling and RNase L inhibition; ideal for translation and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	ABCE1 Rabbit mAb
Remarks/Alias	RLI; OABP; RLI1; ABC38; RNS4I; RNASEL1; RNASELI; ABCE1
Species	Human
Gene ID (Human)	6059
Gene ID	6059
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 500-599 of human ABCE1 (NP_002931.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61221
Protein Weight	67kDa
Shipping	Ice bag