


Anti-Abce1 Rabbit Polyclonal Antibody For Wb And Elisa - P61221

Item Number: CM0025758



Introduction

ABCE1 Rabbit Polyclonal Antibody (P61221) detects human and mouse ABCE1, a ribosome recycling factor and RNase L inhibitor. Suitable for WB and ELISA applications. Ideal for studying translation, mRNA surveillance, and antiviral pathways.

[Learn More](#)

Parameter	Value
Product Name	ABCE1 Rabbit pAb
Remarks/Alias	RLI; OABP; RLI1; ABC38; RNS4I; RNASEL1; RNASELI; ABCE1
Species	Human
Gene ID (Human)	6059
Gene ID	6059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-599 of human ABCE1 (NP_002931.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61221
Protein Weight	67kDa
Shipping	Ice bag