


Anti-Abcf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8Ne71

Item Number: CM0003418



Introduction

High-quality rabbit monoclonal Anti-ABCF1 antibody validated for WB, IHC-P, and ELISA. Targets human ABCF1/ABC50, a translation initiation factor and ATP-binding cassette protein. Ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	ABCF1 Rabbit mAb
Remarks/Alias	ABC27; ABC50; ABCF1
Species	Human
Gene ID	23
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 746-845 of human ABCF1 (Q8NE71).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q8NE71
Protein Weight	96kDa
Shipping	Ice bag