

Anti-Abcf2 Polyclonal Antibody For Wb - Q9Ug63

Item Number: CM0030313

Introduction

High-quality rabbit polyclonal antibody against ABCF2 (SWISS: Q9UG63), validated for WB, IF/ICC, and ELISA. Reacts with human and mouse samples. Ideal for ABC transporter research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	ABCF2 Rabbit pAb
Remarks/Alias	ABC28; HUSSY18; HUSSY-18; EST133090; ABCF2
Species	Human
GeneID (Human)	10061
GeneID	10061
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ABCF2 (NP_005683.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UG63
Protein Weight	71kDa
Shipping	Ice bag