

Anti-Abcg1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P45844

Item Number: CM0030897



Introduction

Rabbit polyclonal antibody targeting human ABCG1, validated for Western blot, immunofluorescence, IHC, and ELISA. Detects ABCG1 in human, mouse, and rat samples, useful for cholesterol/phospholipid efflux studies.

[Learn More](#)

Parameter	Value
Product Name	ABCG1 Rabbit pAb
Remarks/Alias	ABCG1; WHITE1; ABCG1
Species	Human
Gene ID (Human)	9619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 6-120 of human ABCG1 (NP_004906.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45844
Protein Weight	76kDa
Shipping	Ice bag