


Anti-Abcg2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Unq0

Item Number: CM0029398



Introduction

CamelBio offers a rabbit polyclonal antibody targeting human ABCG2 (BCRP), validated for WB, IF/ICC, and ELISA. This 72 kDa multidrug transporter plays a key role in drug resistance and metabolite efflux.

[Learn More](#)

Parameter	Value
Product Name	ABCG2 Rabbit pAb
Remarks/Alias	MRX; MXR; ABCP; BCRP; BMDP; MXR1; ABC15; BCRP1; CD338; GOUT1; MXR-1; CDw338; CDw388; UAQTL1; EST157481; ABCG2
Species	Human
Gene ID (Human)	9429
Gene ID	9429
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human ABCG2 (NP_004818.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNQ0
Protein Weight	72kDa
Shipping	Ice bag