

Anti-Abcg2 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q9Unq0

Item Number: CM0004758



Introduction

High-quality Anti-ABCG2 Rabbit Monoclonal Antibody (Q9UNQ0) validated for WB, IF-P, ELISA in Human, Mouse. Ideal for studying drug resistance, urate transport, and stem cell biology.

[Learn More](#)

Parameter	Value
Product Name	ABCG2 Rabbit mAb
Remarks/Alias	MRX; MXR; ABCP; BCRP; BMDP; MXR1; ABC15; BCRP1; CD338; GOUT1; MXR-1; CDw338; CDw388; UAQTL1; EST157481
Species	Human
GeneID (human)	9429
GeneID	9429
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 101-200 of human ABCG2 (NP_004818.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNQ0
Protein Weight	68kDa/72kDa
Shipping	Ice bag