


Anti-Abcg8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H221

Item Number: CM0019794



Introduction

Anti-ABCG8 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human and Rat. Targets sterolin-2, a key sterol transporter involved in dietary cholesterol absorption and biliary excretion. Suitable for IVD research.

[Learn More](#)

Parameter	Value
Product Name	ABCG8 Rabbit pAb
Remarks/Alias	GBD4; STSL; STSL1; ABCG8
Species	Human
Gene ID (Human)	64241
Gene ID	64241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human ABCG8 (NP_071882.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9H221
Protein Weight	76kDa
Shipping	Ice bag