

Anti-Abf3 Rabbit Polyclonal Antibody (Arabidopsis Thaliana) - Wb, Elisa - Abf3

Item Number: CM0023729

Introduction

Rabbit polyclonal antibody targeting Arabidopsis thaliana ABF3 (abscisic acid responsive elements-binding factor 3). Validated for Western blot and ELISA applications. Recombinant immunogen covers amino acids 260-360. Ideal for plant stress signaling research.

[Learn More](#)

Parameter	Value
Product Name	ABF3 Rabbit pAb
Remarks/Alias	abscisic acid responsive elements-binding factor 3; AtABF3; DC3 PROMOTER-BINDING FACTOR 5; DPBF5; F17I5.200; ABF3
Species	Arabidopsis thaliana
Gene ID	829547
Immunogen	Recombinant fusion protein containing amino acids 260-360 of Arabidopsis thaliana ABF3 (NP_567949.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
Shipping	Ice bag