

Anti-Abhd4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tb40

Item Number: CM0018988

Introduction

Rabbit polyclonal antibody against human ABHD4 (alpha/beta hydrolase domain-containing protein 4), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying N-acylphosphatidylethanolamine lipase and endocannabinoid biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	ABHD4 Rabbit pAb
Remarks/Alias	ABH4; ABHD4
Species	Human
GeneID (Human)	63874
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 223-342 of human ABHD4 (NP_071343.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TB40
Protein Weight	39kDa
Shipping	Ice bag