

Anti-Abhd6 Rabbit Mab For Wb, If-F, Ihc-P, Elisa - Q8R2Y0

Item Number: CM0011305

Introduction

KO validated recombinant Abhd6 rabbit monoclonal antibody for WB, IF-F, IHC-P, ELISA. Detects ABHD6 in human, mouse, rat. Ideal for lipid metabolism and endocannabinoid signaling studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Abhd6 Rabbit mAb
Remarks/Alias	AA673485; AV065425; 0610041D24Rik; ABHD6
Species	Mouse
Gene ID (Human)	57406
Gene ID	66082
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 119-315 of mouse Abhd6 (NP_001317993.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-F, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8R2Y0
Protein Weight	38kDa
Shipping	Ice bag